

METHOD OF COLLECTION

PUS-WOUND



IMPORTANT

IF YOU HAVE DIFFERENT TYPES OF SAMPLES TO TAKE, THE BLOOD SAMPLE SHOULD ALWAYS BE TAKEN LAST.

MATERIAL

To perform this sampling, we provide you with:

- 1 swab.
- 1 tube with a pink cap.



IMPORTANT

The tube cannot be rinsed.
It contains an AMIES liquid medium essential for the analysis.

Before the sample is taken:
Clean the wound by rinsing the area with saline solution.



STEP 1


Write your name, surname and the time of collection on the label of the tube provided.

STEP 2


Wash your hands thoroughly.

STEP 3

Open the sterile packaging and remove the swab. Be careful not to touch the tip of the swab with your fingers.

STEP 4

Apply the smear to the edges of the wound. Take care not to touch the skin.

STEP 5

Transfer the swab into the tube and break the swab stem at the red line.

STEP 6

Seal the tube tightly and place it in the bag provided for this purpose.