



SAMPLING **METHODS**

ACIDIFIED 12-HOUR URINE

IF YOU HAVE DIFFERENT TYPES OF SAMPLE TO TAKE, THE ▲ IMPORTANT BLOOD SAMPLE SHOULD ALWAYS BE TAKEN LAST.

EQUIPMENT

To take this sample, you need:

I clean collection container (not provided).

I jar with a blue cap (provided).

▲ IMPORTANT

The supplied iar may not be rinsed.

It contains the chloridric acid essential for the analysis.

BEFORE TAKING YOUR SAMPLE: Please ask your doctor to find out whether your current treatment (medication, dietary supplement, etc.) should be temporarily suspended.



🕥 LIMS-MBnext | Medical analysis laboratory



Write your name, first name and date of birth on the label provided jar

STEP 2



•

Wash your hands thoroughly.

STEP 3

Urinate <u>a first time in the evening in the</u> <u>toilet</u> (e.g. at 8:00 pm), <u>then start</u> <u>collecting your nightly urine for 12 hours</u> (including urine from 8:00 am) in the clean container you have prepared. Store the urine in the refrigerator for the duration of the collection.



After mixing all the collected urine, open the provided jar and transfer the collected urine to fill it.



Close the jar tightly and place it in the provided bag.